BOTTLE ASSEMBLY

EXPRESSION COLLECTION COMBO KIT (HF)

How to Assemble the Membrane Valve:

Hold the Silicone Membrane and insert the small stem into the top-right opening of the Valve.

Press hard on the Membrane to ensure that the stem is completely inserted and is snug against the Valve opening.

Membrane Valve assembly is complete.

Note: Verify that the membrane is flat. A curled membrane will result in suction loss.


Hold the Breastshield Body in one hand and make sure the bottle-thread side opening is facing up.

Insert the Membrane Valve into the bottom of the Breastshield Body.

Confirm that the Membrane Valve is completely installed.

If desired, place the Silicone Cushion into the Flange of the Breastshield.

Push the center of the Silicone Cushion into the Flange.

Secure the tip of the Cushion around the edge of the Flange.

Ensure that the tip of the Cushion is securely in place.

Note: Your Electric Breast Pump will work with or without the Silicone Cushion. Your comfort depends on your preference.

Insert the Flange into the Breastshield Body.

Place the Diaphragm in the Cup of the Breastshield Body.

Note: Do not reverse the order of steps 6-9. Doing so will result in suction loss.

Ensure that the Diaphragm is securely in place in order to achieve proper suction.

Place the Cap on the Cup and Diaphragm.

Note: Align the hinge cover over the hinge. Press tightly on the cap to seal.

The Electric Diaphragm Kit is assembled.

Twist the Bottle Adapter onto the Breastshield Body.

Twist on the Collection Bottle for a secure fit.

Connect the Tubing to the Tubing connector loaded on the Cap.

Note: Make sure there is not moisture in the Tubing.

Push firmly together for a tight seal.

The Expression Collection Combo Kit is assembled.

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