




BP303 (Standard) / BP401 (Melody) / MP202 (Thermometer)

Easy to use, interchangeable pacifiers for your baby.

The award-winning, patented Detachable Pacifier has three versatile, interchangeable attachments that make caring for your baby a joy. The Standard Module is perfect for daily use. The Melody Module soothes your baby with a sweet lullaby. Additionally, the Thermometer Module accurately detects your baby's temperature. All your baby's pacifier needs in one simple solution.

How to Use:

1. Before first use, wash the Detachable Pacifiers thoroughly with soap and hot water to sterilize.
2. Before each use, check for visible wear or damage.
3. Rinse if necessary before giving to your baby.
4. Use alone or attach to your desired Module by following the "How to Change Modules" instructions below.
5. Standard Module or Melody Module:
 - a. Place pacifier into baby's mouth.
 - b. As your baby suckles the pacifier, the melody module will be activated by the air pressure change inside the nipple.
 - c. Rinse and dry the Detachable Pacifier. Be sure to store the modules in the provided storage case when not in use.
6. Thermometer Module:
 - a. Place pacifier in baby's mouth. Leave in baby's mouth for 6 to 8 minutes to allow the nipple to warm up to give a more accurate reading.
 - b. Press the ON/OFF button to turn the thermometer on. The window will display all symbols  and a beep sound will be heard. This indicates that the unit is working. The last measured temperature will then display for approximately two seconds.
 - c. The flashing 0°C/0°F indicates the thermometer is taking your baby's temperature.
 - d. After three to four minutes, ten beeps will be heard if your baby's temperature is within the normal range. If ten triple-beeps are heard, that indicates an elevated temperature.
 - e. The thermometer will turn off automatically after 10 minutes. However, to prolong battery life, it is recommended to press the ON/OFF button to turn the thermometer off after each use.
 - f. Repeat operation procedures if interrupted during use.
 - g. Rinse and dry the Detachable Pacifier. Be sure to store the Thermometer Module in the provided storage case

when not in use.

How to Change Modules:


1. To UNLOCK

- a. Firmly grasp the Detachable Pacifier in one hand by the shield (not the nipple).
- b. With the other hand, hold the Module (not the handle) and turn counter clock-wise 90° until it stops turning.
- c. Pull apart.
- d. Choose appropriate Interchangeable Module.

2. To LOCK

- a. Align the locking tab to reattach. Fully insert the Module into the Detachable Pacifier. If using the Thermometer Module, gently thread the thermometer probe through the hole first.
- b. Turn the Module (not the handle) clock-wise 90° until you feel it "snap" into place.
- c. If Module is not turning either way, the Module is not inserted completely. Remove Module from Detachable Pacifier and repeat the above process.

Warning/Caution:

- Do not tie pacifier around child's neck, as it presents a strangulation danger.
- Avoid eating or drinking at least 1/2 hour before taking temperature.
- Inspect nipple before and after each use. Dispose of and replace the Detachable Pacifier after any sign of visible wear or damage.
- Do not bend or tamper with the thermometer probe of the Thermometer Module. Only use the Thermometer Module as needed to measure baby's temperature.
- After washing the Detachable Pacifier, be sure to drain the nipple component of the Detachable Pacifier of any water and allow to dry completely before attaching to any of the Interchangeable Modules to protect the electronic components from mildew and other problems. Do not let the Modules come in contact with any chemical solvents.
- Do not expose the Interchangeable Modules to extreme heat or cold. Do not expose the Interchangeable Modules to direct sunlight. Do not get the Interchangeable Modules wet. To clean the Interchangeable Modules, wipe with a damp cloth.
- Store your Thermometer Module in a cool, dry place away from children in the storage box provided. Use only with adult supervision.
- The Thermometer Module is a medical device and is not a toy. It should only be used when taking a temperature. However, we encourage your child to use the Detachable Pacifier regularly as it will help your child accept the pacifier when taking their temperature.
- The Thermometer Module is specifically designed for use with the Detachable Pacifier and should only be used orally. Do not use the Thermometer Module by any other means than what is explained in the instructions. The readings reflect baby's oral temperature which is often 0.5°C/1.0°F lower than actual body temperature.
- A Battery symbol () will display on the Thermometer Module when battery is almost exhausted. This indicates that the Module needs to be replaced. Do not disassemble the unit. Dispose of the entire Module after battery dies, batteries are non-replaceable. Do not dispose of any Module containing a battery in or near a fire as it may explode.
- Components, such as the battery, may be harmful or fatal if swallowed. Immediately consult a doctor.

Specifications:

BPA free.

Medium size cherry style silicone nipple.

Detachable Pacifiers can be sterilized in boiling water.

Display unit: 0.1°C/0.1°F.

Accuracy: ±0.1°C/0.2°F between 35.0°C/95.0°F ~ 38.0°C/100.4°F. Accuracy ±0.2°C/0.4°F for remainder of range at room temperature of 22.0°C/71.6°F.

Measuring range: The Thermometer Module has a Measuring range of 32.0°C/90.0°F ~ 42.9°C/109.9°F.

Outside measuring range will display Lo°C/Lo°F or Hi°C/Hi°F.

Battery: 1.55 VDC (LR41/SR41). Battery life has more than 200 operating hours, approximately 2 years if used for 10 minutes per day. Power consumption 0.2mW.

Module safety locks protect against electronic component damage.

Conforms to the CPSP's (Consumer Product Safety Commission) requirements for pacifiers 16CFR

Part 1511.

Patent:

USA 5,743,647

Rumble Tuff® trademark owned by

A Cute Baby, Inc.

865 N 1430 W Orem, UT 84057

Ph: 1-801-609-8168 Fax: 1-801-796-2688

e-mail: customerservice@acutebaby.com

www.rumbletuff.com www.acutebaby.com

Made In China

75-ST10800-RB01